# babyScale

minilandbaby.com

**miniland**  
BABY

<b>Español</b>	3	<b>Deutsch</b>	16
<b>English</b>	6	<b>Italiano</b>	20
<b>Português</b>	9	<b>Polski</b>	23
<b>Français</b>	13	<b>Русский</b>	27

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# ENGLISH

## TABLE OF CONTENTS

1. INTRODUCTION
2. CONTENT
3. SAFETY INSTRUCTIONS
4. PRODUCT CHARACTERISTICS
5. INSTRUCTIONS FOR USE
6. MAINTENANCE
7. TECHNICAL SPECIFICATIONS
8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

### ■ 1. INTRODUCTION

We would like to congratulate you on your purchase of babyScale, scales for babies and children of up to 20 kg in weight. Read the instructions manual carefully prior to using the apparatus. Keep this manual for future use as it contains important information.

**The characteristics described in this user guide are subject to change with no prior warning.**

### ■ 2. CONTENT

- 1 set of scales
- 3 AAA batteries
- 1 instructions manual and guarantee
- **If any of the aforementioned contents are missing, contact your distributor.**

### ■ 3. SAFETY INSTRUCTIONS

Read the safety instructions carefully prior to using this apparatus and keep it for future reference.

1. You should not strike the scales with sharp or hard objects and avoid strong vibrations.
2. Avoid heavy objects which surpass the maximum capacity of the scales  
20kg/44lb/705.48oz/20000g
3. Never wet the scales. The scales should only be cleaned with a smooth cloth and natural cleaning products.
4. Protect the scales well. If the surface is damaged the readings may be inaccurate.
5. Replace the batteries when the 'LO' indicator appears on the screen.
6. The batteries should be inserted with the poles in the correct position.
7. Remove the batteries before storing the scales for long periods of time.
8. Do not attempt to open or disassemble any part of the device except to replace the batteries.
9. Use the device exclusively for the purpose it is intended for.
10. Do not attempt to repair this product on your own, as opening and handling it may damage the device or expose it to hazards.

**IMPORTANT: Please note that any change or modification to the equipment not performed by the Miniland technical service will not be covered by the product guarantee.**

### ■ 4. PRODUCT CHARACTERISTICS

#### 4.1. DESCRIPTION OF THE SCALES (Figure 1)

- a. The housing
- b. Weight unit button kg/lb/oz/g
- c. "HOLD" button
- d. LCD screen
- e. Tare button

- f. ON/OFF button
- g. Battery compartment
- h. Brackets

#### 4.2. DISPLAY (Figure 2)

1. "HOLD" function on/off icon
2. Battery icon
3. Minus sign
4. Tare icon
5. Digits for the weight measurement
6. Measurement unit icons

### ■ 5. INSTRUCTIONS FOR USE

#### 5.1 PREPARATION

##### 5.1.1. CHANGING THE BATTERIES (Figure 3)

Change the batteries whenever the image on the screen is unclear and the word "LO" appears on the display.

1. Remove the cover to the battery compartment (g) located in the scales.
2. Replace the AAA type batteries. Check the batteries are correctly in place in accordance with the diagram in Figure 3.
3. Close the cover to the compartment correctly.

##### 5.1.2. OPERATING ENVIRONMENT

Be sure to place the scales on a firm horizontal surface in order to avoid damage to the equipment and to ensure the accuracy of the readings.

#### 5.2 USING THE SCALES

1. Press the "ON/OFF" button (f) to turn the scales on. All the display (d) items will appear for a few seconds.
2. Press the "Weight Unit" button (b) when "0.00" appears on the screen to set the weight to "kg", "lb", "oz" or "g". The option you have selected will appear on the display (6).
3. The scales are now ready for use.
4. The Tare function enables you to remove the weight of accessories such as towels to prevent registering a false reading when weighing your baby in this manner. You should first weigh these objects and then press the Tare button (e). The word "TARE" will appear on the display and only what is placed on the scales from then on will be weighed.
5. Place the baby on the scales when "0.00" (kg/lb/oz/g) appears on the screen.
6. The definitive weight will be registered on the display when the word "HOLD" appears. The scales will turn off automatically after one minute.
7. If the baby moves during the process you may use the "HOLD" function. Turn on the scales (conducting the tare process where necessary) and place the baby in position. You will notice the reading on the display varies slightly when the baby moves. Press the "HOLD" button (c) and the word "HOLD" will start to flash on the display. As long as the word "HOLD" is not flashing the scales will read your baby's weight correctly. The scales will turn off automatically after one minute.

**Note:** If you place an object on the scales before turning them on, the reading will be "0.00" when they are turned on (due to the fact the scales conduct an automatic tare when they are turned on). Remove the object and press the "Tare" button (e) to give the empty scales a reading of "0.00".

## ENGLISH

### 5.3. IDENTIFYING AND SOLVING PROBLEMS

PROBLEM	POSSIBLE CAUSES	SOLUTION
The scales will not turn on	The batteries in your scales are flat	Replace the batteries as illustrated in item 5.1.1.
	The batteries have not been put in place correctly	Insert the batteries in accordance with the instructions in item 5.1.1.

## 6. MAINTENANCE

### 6.1. CLEANING

- Clean the surface of the scales with a smooth cloth and natural cleaning products (when necessary).
- Do not use solvents or abrasive detergents.
- Never wet the scales.

## 7. TECHNICAL SPECIFICATIONS

- Units of weight: kg / lb / oz / g
- Maximum capacity: 20 kg / 44 lb / 705.48 oz / 20000 g
- Intervals: 0.01 kg / 0.02 lb / 0.1 oz / 1 g
- Accuracy: 0.6% (of the weight of the object)  $\pm 1d$ , whereby  $d = 10g$

1kg  $\pm 16g$

2kg  $\pm 22g$

3kg  $\pm 28g$

4kg  $\pm 34g$

5kg  $\pm 40g$

6kg  $\pm 46g$

7kg  $\pm 52g$

8kg  $\pm 58g$

9kg  $\pm 64g$

10kg  $\pm 70g$

15kg  $\pm 100g$

20kg  $\pm 130g$

- Overload indicator: "EEEE"
- Low battery indicator: "LO"
- 1-minute automatic off function
- Batteries: 3 x AAA 1.5V

## 8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

- Remember to remove the batteries and dispose of them in a responsible manner when disposing of the baby monitor.
  - You should not mix products bearing the crossed-out rubbish bin symbol with your domestic rubbish.
  - These products should be taken to the collection points allocated by the local authorities for appropriate treatment.
- Alternatively you may contact the point of sale where you purchased the product.

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**We declare under our own responsibility that the product:**

babyScale · Miniland Baby · 89187 · Miniland, S.A.

**to which this declaration refers conforms with the relevant standards or other standardising documents:**

- **Electromagnetic Compatibility (EMC) Directive 2014/30/EU.**
  - EN 55014-1:2006+A1:2009+A2:2011 Electromagnetic compatibility — Requirements for house hold appliances, electric tools and similar apparatus — Part 1: Emission.
  - EN 55014-2:1997+A1:2001+A2:2008 Electromagnetic compatibility — Requirements for house hold appliances, electric tools and similar apparatus — Part 2: Immunity — Product family standard.
- **Restriction of Hazardous Substances (RoHS) Directive 2011/65/EU.**
- **REACH (Registration, Evaluation, Authorization and Restriction of Chemical substances) Regulation (EC) No 1907/2006.**
- **Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU.**

**Miniland, S.A.**  
Chief Executive Officer  
José Verdú Francés  
Onil, Alicante, Spain.  
11th July 2016

Miniland S.A. P.Ind. La Marjal I C/ La Patronal 8-10  
03430 ONIL (Alicante) SPAIN

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Made in China. Keep this information.

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© Miniland, S.A. 2016